

# Spinneret Option


## CLIP SPINNERET

This Standard Spinneret is Easy to Handle.



(Applicable Equipment) NF Series, NANON Series, SD Series



**Clip Spinneret**



(ex. setting)



### Major Features

#### ◆ Easy Handling

- Easily washed  
Only metal connector should be washed. It can be easily washed by putting in the bottle with solvents.
- Easy setting  
For each spinning, a set of nozzle, metal connector, tubing and syringe is required to attach or detach to the clip spinneret. You do not need to detach or setting the clip spinneret so that you can make samples effectively.

The clip spinneret is useful tool to make nanofibers samples from various types of solutions in a shot time.